



Miis



horus SCOPE DDC 200 **DIGITAL DERMOSCOPE**

- Examine dermis/epi-dermis and cool/warm color with one single tool
- Improve quality and efficiency of making the rounds and even reach those in remote areas.
- Minimize paperwork with images or videos taken, stored, and processed digitally.






DIGITAL DERMOSCOPE

HORUS DDC200


POLARIZED

warm color

cool color



Melanoma




Melanoma


NON-POLARIZED

warm color

cool color




Melanoma




Melanoma


clinical showcase




Alopecia Areata



Basal cell carcinoma



Subungual hemorrhage



Blue venus

LENS

Dimension	117x90x203mm
Weight	350g
View Angle	20mm(Diagonal)
Camera/video flash light	natural white LED
Light polarization Function	Change the polarized light function by the switch button on the lens

CONTROL UNIT

Focus	Auto/Manual Focus
Image Resolution	2560x1920 pixels
Video Resolution	1280x960 pixels
LCD Monitor	3.5" Full Color TFT-LCD touch screen
Image Format	JPEG(Photograph)/H.264(Video)
Weight(Including battery)	238g(285g)

RECHARGEABLE

External Power	Power 100~240vac, 50/60Hz
Operating time	3 hours
Charging time	5 hours



Medimaging Integrated Solution Inc.
晉弘科技股份有限公司

• Tel: +886-3-5798860 • E-mail: sales@miis.com.tw • www.miis.com.tw
• 6F, No. 11, Park Avenue II, Hsinchu Science Park, Hsinchu 300091, Taiwan, R.O.C.

